EIRL Viewer DICOM適合性宣言書 / EIRL Viewer DICOM
Conformance Statement
Software Version 2.1, 2.2

ILR-E03001-085 Revision 01

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The supported network DICOM service (SOP) classes are divided into three categories and implemented with their respective roles (User/Provider):

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This document is intended for readers who have a fundamental understanding of the DICOM standard. The target audience includes:

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This document is a conformance statement for the NEMA "Digital Imaging and Communication in Medicine" (DICOM 3.0) standard.

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This document alone does not ensure the interoperability of the EIRL Viewer with other systems.

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Interoperability testing with other systems must be conducted prior to deployment.

3.4. DEFINITIONS, TERMS AND ABBREVIATIONS

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3.5. REFERENCES

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard Available free at <http://medical.nema.org/>

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The EIRL Viewer AE functions as a Storage Service Class SCP, accepting C-STORE-RQ from authorized Storage SCU AEs and receiving various DICOM objects (images, presentation states, key objects), and is capable of registering and storing them within the device.

4.1.3. Sequencing of Real World Activities

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When a user sends images from PACS to the EIRL Viewer AE, the images are stored in EIRL Viewer and become available for viewing.

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The EIRL Viewer AE provides standard conformance for the SOP classes listed in Table 1.

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SOP Class Name	SOP Class UID	SCU	SCP
CT Image Storage SOP Class	1.2.840.10008.5.1.4.1.1.2	No	Yes
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	No	Yes
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	No	Yes
Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	No	Yes
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	No	Yes
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	No	Yes

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Upon device startup, the EIRL Viewer AE automatically initializes as a DICOM communication process and enters a listening state for incoming connections.

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Concurrent associations can be established with multiple Storage SCU AEs. The maximum number of simultaneous associations is specified in Table 3. However, this limit may be constrained by available computer system resources.

Table 3. EIRL Viewer AE の受諾側アソシエーションの数 / Number of Acceptor Associations for EIRL Viewer AE

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The EIRL Viewer AE does not support asynchronous communication.

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The EIRL Viewer AE contains the following implementation identifying information.

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Implementation Class UID	1.3.6.1.4.1.50633.1
Implementation Version Name	X.X.X.X

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X.X.X.X represents the EIRL Viewer software version.

4.2.1.3. Association Initiation Policy

N/A

4.2.1.4. Association Acceptance Policy

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4.2.1.4.1.1. Description and Sequencing of Activities

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The EIRL Viewer AE accepts associations exclusively from Storage SCUs.

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Image display can be performed by transmitting CT/MR/CR/DX/SC images to the EIRL Viewer AE.

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名前 / Name	UID
JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1]): Default Transfer Syntax for Lossless JPEG Image Compression	1.2.840.10008.1.2.4.70
JPEG Baseline (Process 1): Default Transfer Syntax for Lossy JPEG 8 Bit Image Compression	1.2.840.10008.1.2.4.50
JPEG Extended (Process 2 & 4): Default Transfer Syntax for Lossy JPEG 12 Bit Image Compression (Process 4 only)	1.2.840.10008.1.2.4.51
Explicit VR Little Endian	1.2.840.10008.1.2.1
Implicit VR Little Endian	1.2.840.10008.1.2
Explicit VR Big Endian	1.2.840.10008.1.2.2

4.3. NETWORK INTERFACES

4.3.1. Supported Communications Stacks (parts 8, 9)

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The following extended character set is supported:

- ISO_IR 100

Chapter 7. SECURITY

N/A

Chapter 8. ANNEXES

8.1. IOD Contents

8.1.1. Created SOP Instance(s)

N/A

8.1.2. Usage of Attributes From Received IODs

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The following attributes are utilized when EIRL Viewer displays received images.

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