EIRL DICOM適合性宣言書 / EIRL DICOM Conformance Statement

ILR-ELB001-001 Revision 12

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プログラム	プログラムのバージョン
Program	Program Version
EIRL Brain Aneurysm	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Brain Metry	1.13, 2.0以降 / 1.13 and 2.0 or later
EIRL Brain Segmentation	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest Nodule	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest XR	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest Metry	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest CT Basic	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest CT	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later
EIRL Chest TB	1.12, 1.13, 2.0以降 / 1.12, 1.13 and 2.0 or later

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The supported network DICOM service (SOP) classes are organized into the following three categories, with specific roles (User/Provider):

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Services can be specified as SCU (Service Class User), SCP (Service Class Provider), or as optional.

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The UID values for each SOP class are classified as follows:

UID Value	UID Name	カテゴリ Category
1.2.840.10008.4.2	Storage Service Class	転送 / Transfer
1.2.840.10008.5.1.4.1.1.4	MR Image Storage	転送 / Transfer
1.2.840.10008.5.1.4.1.1.2	CT Image Storage	転送 / Transfer
1.2.840.10008.5.1.4.1.1.1	CR Image Storage	転送 / Transfer
1.2.840.10008.5.1.4.1.1.1.1	Digital X-Ray Image Storage - For Presentation	転送 / Transfer
1.2.840.10008.5.1.4.1.1.88.67	X-Ray Radiation Dose SR Storage	転送 / Transfer
1.2.840.10008.5.1.4.1.1.11.1	Grayscale Softcopy Presentation State Storage SOP Class	転送 / Transfer
1.2.840.10008.5.1.4.1.1.7	Secondary Capture Image Storage	転送 / Transfer

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The supported media storage application profiles are categorized into three types with their respective roles.

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N/A

Chapter 2. TABLE OF CONTENTS

1. CONFORMANCE STATEMENT OVERVIEW
2. TABLE OF CONTENTS
3. INTRODUCTION
 - 3.1. REVISION HISTORY
 - 3.2. AUDIENCE
 - 3.3. REMARKS
 - 3.4. DEFINITIONS, TERMS AND ABBREVIATIONS
 - 3.5. REFERENCES
4. NETWORKING
 - 4.1. IMPLEMENTATION MODEL
 - 4.1.1. Application Data Flow Diagram
 - 4.1.2. Functional definitions of AE's
 - 4.1.3. Sequencing of Real World Activities
 - 4.2. AE SPECIFICATIONS
 - 4.2.1. "EIRL" - Specification
 - 4.3. NETWORK INTERFACES
 - 4.3.1. Supported Communications Stacks (parts 8, 9)
 - 4.3.2. TCP/IP Stacks
 - 4.4. CONFIGURATION
5. MEDIA INTERCHANGE
 - 5.1. IMPLEMENTATION MODEL
 - 5.1.1. Application Data Flow Diagram
 - 5.1.2. Functional definitions of AE's
 - 5.1.3. Sequencing of Real World Activities
 - 5.1.4. File Meta Information for Implementation Class and Version
 - 5.2. AE SPECIFICATIONS
 - 5.2.1. "EIRL" - Specification
 - 5.3. AUGMENTED AND PRIVATE APPLICATION PROFILES
 - 5.4. MEDIA CONFIGURATION
6. SUPPORT OF CHARACTER SETS
7. SECURITY
8. ANNEXES
 - 8.1. IOD Contents
 - 8.1.1. Created SOP Instance(s)
 - 8.1.2. Usage of Attributes From Received IODs
 - 8.2. Data Dictionary of Private Attributes
 - 8.2.1 Private Group LPIXEL 003

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This document is intended for readers who have a basic knowledge of the DICOM standard.

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This document is a conformance statement for NEMA's "Digital Imaging and Communications in Medicine" (DICOM 3.0) standard.

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This document alone does not guarantee interoperability between EIRL and other systems.

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Please ensure interoperability testing with other systems is performed before deployment.

3.4. DEFINITIONS, TERMS AND ABBREVIATIONS

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3.5. REFERENCES

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard

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The EIRL AE waits in the background to receive connection requests. When an association with Storage Service Class SOP classes is established, it stores the received instances in local storage.

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Next, from the received studies, it performs the following image analyses on series that meet the applicability conditions for each application:

- 画像解析（脳動脈瘤読影支援） / Image analysis (brain aneurysm detection support)
- 脳計測（Evans Index） / Brain measurement (Evans Index)
- 脳計測（Callosal Angle） / Brain measurement (Callosal Angle)
- 脳計測（白質高信号領域） / Brain measurement (white matter hyperintensity regions)
- 画像解析（胸部X線肺結節読影支援） / Image analysis (chest X-ray lung nodule detection support)
- 関心領域計測（肺計測） / Region of interest measurement (lung measurement)
- 関心領域抽出（頭部CTセグメンテーション） / Region of interest extraction (head CT segmentation)
- 関心領域抽出（胸部CTセグメンテーション） / Region of interest extraction (chest CT segmentation)
- 肺結節(候補)領域検出（Chest CT2） / Lung nodule (candidate) region detection (Chest CT2)
- 画像解析（Chest XR） / Image analysis (Chest XR)
- 画像解析（Chest TB） / Image analysis (Chest TB)

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Finally, after image analysis, it transmits the analysis results as DICOM GSPS/SC/Overlay:

- 検出された候補点の座標（脳動脈瘤読影支援、胸部X線肺結節読影支援、ChestXR、ChestTB）
Coordinates of detected candidate points (brain aneurysm detection support, chest X-ray lung nodule detection support, ChestXR, ChestTB)
- 過去比較機能で、過去に検出されたが検出されなくなった点の座標（胸部X線肺結節読影支援）
Coordinates of points that were previously detected but are no longer detected in temporal comparison (chest X-ray lung nodule detection support)
- 最大側脳室前角幅とその位置での大脳幅を表す線 (Evans Index)
Lines representing the maximum frontal horn width and cerebral width at that position (Evans Index)
- 側脳室の上内面を線で近似する線 (Callosal Angle)
Lines approximating the superior medial surface of the lateral ventricles (Callosal Angle)
- 白質高信号領域を塗りつぶす
Fill white matter hyperintensity regions
- 計測された関心領域の輪郭線を描画 (肺計測)
Draw contour lines of measured regions of interest (lung measurement)
- 計測された関心領域を塗りつぶす（頭部CTセグメンテーション）
Fill measured regions of interest (head CT segmentation)
- 計測された関心領域を塗りつぶす（胸部CTセグメンテーション）
Fill measured regions of interest (chest CT segmentation)

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EIRL can save analysis results in any of the following formats:

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Additionally, automatic generation of MIP images is available for brain aneurysm detection support. In that case, the generated MIP images are saved as SC and GSPS.

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Additionally, automatic generation of RaySum images (coronal / sagittal) is available for Chest CT segmentation/Chest CT2. In that case, the generated RaySum images are saved as SC and GSPS.

4.1.3. Sequencing of Real World Activities

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No specific sequencing requirements for real world activities.

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The EIRL AE provides standard conformance to the SOP classes listed in [Table 1](#).

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SOP Class Name	SOP Class UID	SCU	SCP
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Yes	No

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When this device is started, the EIRL AE automatically launches as a DICOM communication process and enters a receiving standby state.

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Associations can be established simultaneously with multiple Storage SCU AEs. The maximum number of simultaneous associations is shown in Table 3. However, this maximum may be limited by computer system resources.

Table 3. EIRL AE の受諾側アソシエーションの数 / Number of Acceptor Associations for EIRL AE

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The EIRL AE does not provide asynchronous communication.

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The EIRL AE has the following implementation identifying information.

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Implementation Class UID	1.3.6.1.4.1.50633.1
Software Versions	X.X.X.X

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4.2.1.3. Association Initiation Policy

4.2.1.3.1. Activity: GSPS/SC/Overlayの保存 / Storage of GSPS/SC/Overlay

4.2.1.3.1.1. Description and Sequencing of Activities

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Required Pre-registration Information for Storage SCP AE

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After registering the above information, EIRL analysis results are transferred to the registered Storage SCP.

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名前	UID
Explicit VR Little Endian	1.2.840.10008.1.2.1
Implicit VR Little Endian	1.2.840.10008.1.2
Explicit VR Big Endian	1.2.840.10008.1.2.2

4.2.1.4. Association Acceptance Policy

4.2.1.4.1. Activity: 画像の解析 / Image Analysis

4.2.1.4.1.1. Description and Sequencing of Activities

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The EIRL AE accepts associations only from Storage SCUs.

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Image analysis can be performed by sending MR/CR/DX images to the EIRL AE.

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名前	UID
JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1]): Default Transfer Syntax for Lossless JPEG Image Compression	1.2.840.10008.1.2.4.70
JPEG Baseline (Process 1): Default Transfer Syntax for Lossy JPEG 8 Bit Image Compression	1.2.840.10008.1.2.4.50
JPEG Extended (Process 2 & 4): Default Transfer Syntax for Lossy JPEG 12 Bit Image Compression (Process 4 only)	1.2.840.10008.1.2.4.51
JPEG 2000 Image Compression (Lossless Only)	1.2.840.10008.1.2.4.90
Explicit VR Little Endian	1.2.840.10008.1.2.1
Implicit VR Little Endian	1.2.840.10008.1.2
Explicit VR Big Endian	1.2.840.10008.1.2.2

4.3. NETWORK INTERFACES

4.3.1. Supported Communications Stacks (parts 8, 9)

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Provides the DICOM upper layer protocol for TCP/IP as defined in DICOM PS3.8.

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The integration BOX supports the following network interface configurations.

Ethernet 1000Base-TX is the recommended configuration, with 100Base-TX available as an option.

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The EIRL AE can perform the following functions:

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5.1.4. File Meta Information for Implementation Class and Version

File Meta Information Version	00 01
Implementation Class UID	1.3.6.1.4.1.50633.1
Software Versions	X.X.X.X

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EIRL provides standard conformance to the DICOM Interchange Option of the Media Storage Service Class.

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The supported application profiles and roles are shown in the table below.

Application Profiles Supported	Real World Activity	Role	SC Option
STD-GEN-CD	Load directory or file	FSR FSC	Interchange

5.3. AUGMENTED AND PRIVATE APPLICATION PROFILES

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The following extended character sets are supported.

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The following table uses abbreviations. These abbreviations are used in the "Presence of Module" column.

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The following abbreviations are used in the "Source" column.

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Modules for Generated Grayscale Softcopy Presentation State SOP Instance IOD:

IE	モジュール / Module	モジュールの存在 / Presence of Module
Patient	Patient	ALWAYS
Study	General Study	ALWAYS
	Patient Study	ANAP
Series	General Series	ALWAYS
	Presentation Series	ALWAYS
Equipment	General Equipment	ALWAYS
Presentation State	Presentation State Identification	ALWAYS
	Presentation State Relationship	ALWAYS
	Presentation State Shutter	ALWAYS
	Presentation State Mask	ALWAYS
	Displayed Area	ALWAYS
	Graphic Annotation	ALWAYS
	Graphic Layer	ALWAYS
	Softcopy Presentation LUT	ALWAYS
	SOP Common	ALWAYS

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Presentation State Relationship Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
Referenced Series Sequence	(0008,1115)	SQ	元画像のUIDを参照 Reference to original image UID	ALWAYS	SRC

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Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Graphic Annotation Units	(0070,0005)	CS	固定値: Fixed value: "PIXEL"	ALWAYS	AUTO
>> Graphic Dimensions	(0070,0020)	US	固定値: Fixed value: "2"	ALWAYS	AUTO
>> Number of Graphic Points	(0070,0021)	US	固定値: Fixed value: "2"	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Graphic Data	(0070,0022)	FL	候補点のx,y座標と、円形の右側の点のx,y座標 (脳動脈瘤) x,y coordinates of candidate points and right-side points of circles (brain aneurysm) 又は or 最大側脳室前角幅とその位置での大脳幅を表す線の起点と終点 (Evans Index) Start and end points of lines representing maximum ventricular frontal horn width and cerebral width at that position (Evans Index) 又は or 側脳室の上内面を線で近似する線の起点と終点 (Callosal Angle) Start and end points of lines approximating the superior medial surface of lateral ventricles (Callosal Angle) 又は or 白質高信号領域を囲む多角形の点のx,y座標 x,y coordinates of polygon points surrounding white matter hyperintensity regions 又は or 検出領域を示す長方形の左上及び右下のx,y座標 (肺結節) x,y coordinates of top-left and bottom-right corners of rectangles indicating detection regions (pulmonary nodules) 又は or 過去比較機能で、過去に検出されたが検出されなくなった点を示すx,y座標と、円形の右側の点のx,y座標 (肺結節) x,y coordinates indicating points detected in the past but no longer detected by the past comparison function, and x,y coordinates	ALWAYS	ANALYSIS

Attribute Name	Tag	VR	Value	Presence of Value	Source
同上 Same as above	同上 Same as above	同上 Same as above	つづき Continuation 頭蓋骨内で周囲より高吸収及び低吸収の領域を示す多角形の点のx,y座標 (頭部CTセグメンテーション) x,y coordinates of polygon points indicating high-attenuation and low-attenuation regions within the skull (head CT segmentation) 又は or 関心領域を示す多角形の点のx,y座標 (胸部CTセグメンテーション) x,y coordinates of polygon points indicating regions of interest (chest CT segmentation) 又は or 肺結節(候補)領域を示す多角形の点のx,y座標 (Chest CT2) x,y coordinates of polygon points indicating pulmonary nodule (candidate) regions (Chest CT2) 又は or 2D 長径・短径を示す測定線分の起点及び終点のx,y座標 (Chest CT segmentation/Chest CT2) x,y coordinates of the start and end points of line segments indicating 2D long diameter (LAD) / short diameter (SAD) measurements (Chest CT segmentation/Chest CT2) 又は or 各領域に対応する塗りつぶし円の中心点のx,y座標と、円周上の1点のx,y座標 (Chest CT segmentation/Chest CT2 RaySum) x,y coordinates of the center point and one point on the circumference of filled circles	同上 Same as above	同上 Same as above

Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Graphic Type	(0070,0023)	CS	固定値: Fixed value: "CIRCLE" (脳動脈瘤、肺結節、肺計測、Chest CT segmentation/Chest CT2 RaySum) Fixed value: "CIRCLE" (brain aneurysm, pulmonary nodules, lung measurement, Chest CT segmentation/Chest CT2 RaySum) 又は or "POLYLINE" (脳計測、白質高信号領域、肺結節、肺計測、頭部CTセグメンテーション、胸部CTセグメンテーション、Chest CT2、Chest XR、Chest TB) "POLYLINE" (brain measurement, white matter hyperintensity regions, pulmonary nodules, lung measurement, head CT segmentation, chest CT segmentation, Chest CT2, Chest XR, Chest TB)	ALWAYS	AUTO
>> Graphic Filled	(0070,0024)	CS	"N" (デフォルト) 又は "Y" (Chest CT segmentation/Chest CT2 RaySum) "N" (default) or "Y" (Chest CT segmentation/Chest CT2 RaySum)	ALWAYS	AUTO
>> Graphic Group ID	(0070,0295)	UL	グループの番号 Group number	ALWAYS	AUTO
> Compound Graphic Sequence	(0070,0209)	SQ	設定でCompound Graphic Sequenceを保存するように指定可能。 Can be configured to save Compound Graphic Sequence in settings.	VNAP	ANALYSIS / CONFIG
>> Compound Graphic Instance ID	(0070,0226)	UL	コンパウンドグラフィックの番号 Compound graphic number	ALWAYS	AUTO
>> Compound Graphic Units	(0070,0282)	CS	固定値: Fixed value: "PIXEL"	ALWAYS	AUTO
>> Graphic Dimensions	(0070,0020)	US	固定値: Fixed value: "2"	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Number of Graphic Points	(0070,0021)	US	固定値: Fixed value: "2"	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Graphic Data	(0070,0022)	FL	候補点を囲む楕円形の左上隅の座標と、右下隅の座標 (脳動脈瘤) Coordinates of top-left and bottom-right corners of ellipse surrounding candidate points (brain aneurysm) 又は or 最大側脳室前角幅とその位置での大脳幅を表す線 (Evans Index) Lines representing maximum ventricular frontal horn width and cerebral width at that position (Evans Index) 又は or 側脳室の上内面を線で近似する線 (Callosal Angle) Lines approximating the superior medial surface of lateral ventricles (Callosal Angle) 又は or 白質高信号領域を囲む多角形の点のx,y座標 x,y coordinates of polygon points surrounding white matter hyperintensity regions 又は or 検出領域を示す長方形の左上及び右下のx,y座標 (肺結節) x,y coordinates of top-left and bottom-right corners of rectangles indicating detection regions (pulmonary nodules) 又は or 過去比較機能で、過去に検出されたが検出されなくなった点を囲む楕円形の左上隅の座標と、右下隅の座標 (肺結節) Coordinates of top-left and bottom-right corners of ellipse surrounding points detected in the past but	ALWAYS	ANALYSIS

Attribute Name	Tag	VR	Value	Presence of Value	Source
同上 Same as above	同上 Same as above	同上 Same as above	つづき Continuation 関心領域を示す多角形の点のx,y座標 (胸部CTセグメンテーション) x,y coordinates of polygon points indicating regions of interest (chest CT segmentation) 又は or 肺結節(候補)領域を示す多角形の点のx,y座標 (Chest CT2) x,y coordinates of polygon points indicating pulmonary nodule (candidate) regions (Chest CT2) 又は or 2D 長径・短径を示す測定線分の起点及び終点のx,y座標 (Chest CT segmentation/Chest CT2) x,y coordinates of the start and end points of line segments indicating 2D long diameter (LAD) / short diameter (SAD) measurements (Chest CT segmentation/Chest CT2) 又は or RaySum投影画像上の各領域に対応する塗りつぶし円を囲む楕円形の左上隅の座標と、右下隅の座標 (Chest CT segmentation/Chest CT2 RaySum) Coordinates of top-left and bottom-right corners of ellipses surrounding filled circles corresponding to each region on RaySum projection images (Chest CT segmentation/Chest CT2 RaySum) 又は or 検出領域を示す長方形の左上及び右下のx,y座標 (Chest XR) x,y coordinates of top-left and bottom-right corners	同上 Same as above	同上 Same as above

Attribute Name	Tag	VR	Value	Presence of Value	Source
>> Compound Graphic Type	(0070,0294)	CS	固定値: Fixed value: "ELLIPSE" (脳動脈瘤、肺結節、肺計測、Chest CT segmentation/Chest CT2 RaySum) Fixed value: "ELLIPSE" (brain aneurysm, pulmonary nodules, lung measurement, Chest CT segmentation/Chest CT2 RaySum) 又は or "MULTILINE" (脳計測、白質高信号領域、肺結節、肺計測、頭部CTセグメンテーション、胸部CTセグメンテーション、Chest CT2、Chest XR、Chest TB) "MULTILINE" (brain measurement, white matter hyperintensity regions, pulmonary nodules, lung measurement, head CT segmentation, chest CT segmentation, Chest CT2, Chest XR, Chest TB)	ALWAYS	AUTO
>> Graphic Filled	(0070,0024)	CS	"N" (デフォルト) 又は "Y" (Chest CT segmentation/Chest CT2 RaySum) "N" (default) or "Y" (Chest CT segmentation/Chest CT2 RaySum)	ALWAYS	AUTO
>> Graphic Group ID	(0070,0295)	UL	グループの番号 Group number	ALWAYS	AUTO

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Graphic Layer Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
Graphic Layer Sequence	(0070,0060)	SQ		ALWAYS	AUTO
> Graphic Layer	(0070,0002)	CS	"LAYER1"	ALWAYS	AUTO
> Graphic Layer Order	(0070,0062)	IS	固定値: Fixed value: "1"	ALWAYS	AUTO
> Graphic Layer Recommended Display Grayscale Value	(0070,0066)	US	固定値: Fixed value: 白 Fixed value: White	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Presence of Value	Source
> Graphic Layer Recommended Display CIELab Value	(0070,0401)	US	固定値: 赤 (デフォルト) Fixed value: Red (default) * Chest CT segmentation (バージョン2.4以降) / Chest CT2 (バージョン2.2以降) は結節サイズに応じた色 (緑 / 黄 / 赤) を各グラフィック・テキストオブ ジェクト単位で適用する Note: Chest CT segmentation (version 2.4 or later) / Chest CT2 (version 2.2 or later) applies size-based per-object colors (green / yellow / red) to individual graphic and text objects	ALWAYS	AUTO
> Graphic Layer Description	(0070,0068)	LO	固定値: Fixed value: "EIRL analysis result"	ALWAYS	AUTO

Softcopy Presentation LUTモジュール:

Softcopy Presentation LUT Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
Presentation LUT Shape	(2050,0020)	CS	固定値: Fixed value: "IDENTITY"	ALWAYS	AUTO

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SOP Common Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
SOP Class UID	(0008,0016)	UI	Grayscale Softcopy Presentation State	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	生成された値 Generated value	ALWAYS	AUTO
Specific Character Set	(0008,0005)	CS	元画像より From original image	ANAP	SRC

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Modules for Generated Secondary Capture Instance IOD:

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IE	モジュール / Module	モジュールの存在 / Presence of Module
Study	General Study	ALWAYS
	Patient Study	ANAP
Series	General Series	ALWAYS
Frame of Reference	Frame of Reference	ANAP
Equipment	General Equipment	ALWAYS
	SC Equipment	ALWAYS
Image	General Image	ALWAYS
	Image Plane	ANAP
	Image Pixel	ALWAYS
	SC Image	ALWAYS
	Modality LUT	ALWAYS
	VOI LUT	ALWAYS
	SOP Common	ALWAYS

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Patient Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
Patient's Name	(0010,0010)	PN	元画像より From original image	VNAP	SRC
Patient ID	(0010,0020)	LO	元画像より From original image	VNAP	SRC
Patient's Birth Date	(0010,0030)	DA	元画像より From original image	VNAP	SRC
Patient Sex	(0010,0040)	CS	元画像より From original image	VNAP	SRC
Referenced Patient Sequence	(0008,1120)	SQ	元画像より From original image	ANAP	SRC
Issuer of Patient ID	(0010,0021)	LO	元画像より From original image	ANAP	SRC
Type of Patient ID	(0010,0022)	CS	元画像より From original image	ANAP	SRC
Issuer of Patient ID Qualifiers Sequence	(0010,0024)	SQ	元画像より From original image	ANAP	SRC
Source Patient Group Identification Sequence	(0010,0026)	SQ	元画像より From original image	ANAP	SRC

Attribute Name	Tag	VR	Value	Presence of Value	Source
Group of Patients Identification Sequence	(0010,0027)	SQ	元画像より From original image	ANAP	SRC
Patient's Birth Time	(0010,0032)	TM	元画像より From original image	ANAP	SRC
Patient's Birth Date in Alternative Calendar	(0010,0033)	LO	元画像より From original image	ANAP	SRC
Patient's Death Date in Alternative Calendar	(0010,0034)	LO	元画像より From original image	ANAP	SRC
Patient's Alternative Calendar	(0010,0035)	CS	元画像より From original image	ANAP	SRC
Quality Control Subject	(0010,0200)	CS	元画像より From original image	ANAP	SRC
Strain Description	(0010,0212)	UC	元画像より From original image	ANAP	SRC
Strain Nomenclature	(0010,0213)	LO	元画像より From original image	ANAP	SRC
Strain Stock Sequence	(0010,0216)	SQ	元画像より From original image	ANAP	SRC
Strain Additional Information	(0010,0218)	UT	元画像より From original image	ANAP	SRC
Strain Code Sequence	(0010,0219)	SQ	元画像より From original image	ANAP	SRC
Genetic Modifications Sequence	(0010,0221)	SQ	元画像より From original image	ANAP	SRC
Other Patient Names	(0010,1001)	PN	元画像より From original image	ANAP	SRC
Other Patient IDs Sequence	(0010,1002)	SQ	元画像より From original image	ANAP	SRC
Referenced Patient Photo Sequence	(0010,1100)	SQ	元画像より From original image	ANAP	SRC
Ethnic Group	(0010,2160)	SH	元画像より From original image	ANAP	SRC
Patient Species Description	(0010,2201)	LO	元画像より From original image	ANAP	SRC
Patient Species Code Sequence	(0010,2202)	SQ	元画像より From original image	ANAP	SRC
Patient Breed Description	(0010,2292)	LO	元画像より From original image	ANAP	SRC
Patient Breed Code Sequence	(0010,2293)	SQ	元画像より From original image	ANAP	SRC

Attribute Name	Tag	VR	Value	Presence of Value	Source
Breed Registration Sequence	(0010,2294)	SQ	元画像より From original image	ANAP	SRC
Responsible Person	(0010,2297)	PN	元画像より From original image	ANAP	SRC
Responsible Person Role	(0010,2298)	CS	元画像より From original image	ANAP	SRC
Responsible Organization	(0010,2299)	LO	元画像より From original image	ANAP	SRC
Patient Comments	(0010,4000)	LT	元画像より From original image	ANAP	SRC
Patient Identity Removed	(0012,0062)	CS	元画像より From original image	ANAP	SRC
De-identification Method	(0012,0063)	LO	元画像より From original image	ANAP	SRC
De-identification Method Code Sequence	(0012,0064)	SQ	元画像より From original image	ANAP	SRC

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General Study Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
Study Instance UID	(0020,000D)	UI	元画像より From original image	ALWAYS	SRC
Study Date	(0008,0020)	DA	元画像より From original image	VNAP	SRC
Study Time	(0008,0030)	TM	元画像より From original image	VNAP	SRC
Referring Physician's Name	(0008,0090)	PN	元画像より From original image	VNAP	SRC
Study ID	(0020,0010)	SH	元画像より From original image	VNAP	SRC
Accession Number	(0008,0050)	SH	元画像より From original image	VNAP	SRC
Study Description	(0008,1030)	LO	元画像より From original image	VNAP	SRC
Issuer of Accession Number Sequence	(0008,0051)	SQ	元画像より From original image	ANAP	SRC
Referring Physician Identification Sequence	(0008,0096)	SQ	元画像より From original image	ANAP	SRC
Consulting Physician's Name	(0008,009C)	PN	元画像より From original image	ANAP	SRC

Attribute Name	Tag	VR	Value	Presence of Value	Source
Consulting Physician Identification Sequence	(0008,009D)	SQ	元画像より From original image	ANAP	SRC
Procedure Code Sequence	(0008,1032)	SQ	元画像より From original image	ANAP	SRC
Physician(s) of Record	(0008,1048)	PN	元画像より From original image	ANAP	SRC
Physician(s) of Record Identification Sequence	(0008,1049)	SQ	元画像より From original image	ANAP	SRC
Name of Physician(s) Reading Study	(0008,1060)	PN	元画像より From original image	ANAP	SRC
Physician(s) Reading Study Identification Sequence	(0008,1062)	SQ	元画像より From original image	ANAP	SRC
Referenced Study Sequence	(0008,1110)	SQ	元画像より From original image	ANAP	SRC
Requesting Service	(0032,1033)	LO	元画像より From original image	ANAP	SRC
Requesting Service Code Sequence	(0032,1034)	SQ	元画像より From original image	ANAP	SRC
Reason For Performed Procedure Code Sequence	(0040,1012)	SQ	元画像より From original image	ANAP	SRC

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SOP Common Module:

Attribute Name	Tag	VR	Value	Presence of Value	Source
SOP Class UID	(0008,0016)	UI	Secondary Capture Image Storage	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	生成された値 Generated value	ALWAYS	AUTO
Specific Character Set	(0008,0005)	CS	元画像より From original image	ANAP	SRC

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For Overlay output, the following IOD modules are added to the original images:

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The overlay output instance is a copy of the original image with the following overlay attributes added, and the SOP Class UID (0008,0016) keeps the value of the original image.

Attribute Name	Tag	VR	Value	Presence of Value	Source
Overlay Rows	(6000,0010)	US	表示エリアの行数 Number of rows in display area	ALWAYS	AUTO
Overlay Columns	(6000,0011)	US	表示エリアの列数 Number of columns in display area	ALWAYS	AUTO
Overlay Type	(6000,0040)	CS	固定値: Fixed value: "G"	ALWAYS	AUTO
Overlay Origin	(6000,0050)	SS	表示エリアの左上の点 Top-left point of display area	ALWAYS	AUTO
Overlay Bits Allocated	(6000,0100)	US	固定値: Fixed value: 1	ALWAYS	AUTO
Overlay Bits Position	(6000,0102)	US	固定値: Fixed value: 0	ALWAYS	AUTO
Overlay Data	(6000,3000)	OW	解析結果を示すビットマップの ピクセルデータ Pixel data of bitmap showing analysis results	ALWAYS	AUTO

8.1.2. Usage of Attributes From Received IODs

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When EIRL analyzes received images, the following attributes are used:

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"LPIXEL 003" provides additional information for the Image Information Entity.

属性名 / Attribute Name	タグ / Tag	VR	VM	Presence of Value	Condition	注記 / Notes
Private Creator	(0021,00xx)	LO	1	ANAP	(*1)	"LPIXEL 003"
Number of Detections	(0021,xx01)	IS	1	ANAP	(*2)(*3)	検出数 Number of detections
Pleural Space Area	(0021,xx10)	IS	1	ANAP	(*2)(*4)	胸腔内空気含有領域の面積 (単位: mm ²) Pleural space area (unit: mm ²)
CP Angle 0	(0021,xx11)	DS	1	ANAP	(*2)(*4)	肋骨横隔膜角/CPアングル0 (単位: °) Costophrenic angle 0 (unit: °)
CP Angle 1	(0021,xx12)	DS	1	ANAP	(*2)(*4)	肋骨横隔膜角/CPアングル1 (単位: °) Costophrenic angle 1 (unit: °)
Cardiothoracic Ratio	(0021,xx13)	DS	1	ANAP	(*2)(*4)	心胸郭比 (単位: %) Cardiothoracic ratio (unit: %)
Heart Width	(0021,xx14)	DS	1	ANAP	(*2)(*4)	心臓幅 (単位: mm) Heart width (unit: mm)
Mediastinum Width Upper	(0021,xx15)	DS	1	ANAP	(*2)(*4)	縦隔幅/上部 (単位: mm) Mediastinum width upper (unit: mm)
Mediastinum Width Middle	(0021,xx16)	DS	1	ANAP	(*2)(*4)	縦隔幅/中部 (単位: mm) Mediastinum width middle (unit: mm)
Mediastinum Width Lower	(0021,xx17)	DS	1	ANAP	(*2)(*4)	縦隔幅/下部 (単位: mm) Mediastinum width lower (unit: mm)
Aortic Arch Diameter	(0021,xx18)	DS	1	ANAP	(*2)(*4)	大動脈弓の径 (単位: mm) Aortic arch diameter (unit: mm)
Longest Axial Diameter	(0021,xx20)	DS	1	ANAP	(*2)(*5)	平均径が最大の領域の2次元長径 (単位: mm) Longest axial diameter of the region with the largest average diameter (unit: mm)
Short Axial Diameter	(0021,xx21)	DS	1	ANAP	(*2)(*5)	平均径が最大の領域の2次元短径 (単位: mm) Short axial diameter of the region with the largest average diameter (unit: mm)

属性名 / Attribute Name	タグ / Tag	VR	VM	Presence of Value	Condition	注記 / Notes
Average Axial Diameter	(0021,xx22)	DS	1	ANAP	(*2)(*5)	平均径が最大の領域の2次元平均径 (単位: mm) Average axial diameter of the region with the largest average diameter (unit: mm)
3D Longest Diameter	(0021,xx23)	DS	1	ANAP	(*2)(*5)	平均径が最大の領域の3次元長径 (単位: mm) 3D longest diameter of the region with the largest average diameter (unit: mm)
Volume	(0021,xx24)	DS	1	ANAP	(*2)(*5)	平均径が最大の領域の体積 (単位: mm ³) Volume of the region with the largest average diameter (unit: mm ³)

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This tag is output only when one or more tags from (0021,xx01 to xx24) are output

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Depending on settings, elements may not be output

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Applicable only for: Brain aneurysm, X-Ray Lung nodule, X-Ray Lung nodule + Chest Metry, Chest XR, Chest XR + Chest Metry, or Chest CT segmentation/Chest CT2

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Applicable only for: Chest Metry, X-Ray Lung nodule + Chest Metry, or Chest XR + Chest Metry. Note that Chest Metry refers to version 2.2 or later.

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Applicable only for: Chest CT segmentation/Chest CT2. Output only to Axial GSPS/SC when one or more regions exist. Note that Chest CT segmentation refers to version 2.4 or later, and Chest CT2 refers to version 2.2 or later.